

Measuring Electron Temperature using Spectral Data and a Collisional Radiative Model for PFRC-II

Christopher Jakuback

August 2020

SULI Summer
Program



PFRC-II Pulse


Thank you:
Dr. Sangeeta Vinoth
Dr. Samuel Cohen

This work was made possible by funding from the Department of Energy for the Summer Undergraduate Laboratory Internship (SULI) program. This work is supported by the US DOE Contract No. DE-AC02-09CH11466.


- **What is PFRC-II?**
- **PFRC-II and the need for Spectral Diagnostics**
- **Use of a Collisional Radiative Model to output Electron Temperature Data**
- **Current Observations**
- **Future Work**



- Princeton Field Reversed Configuration (PFRC)
- “High Beta” plasma containment device
- Relies on an Odd-Parity Rotating Magnetic Field
- Comparatively easy to maintain and construct
- Potentially viable for advanced fusion fuel cycles.



- Visible light from PFRC-II can be used to analyze certain parameters in the plasma.
- Relating net emissions of the Balmer Series provides insight into the plasma over the course of a pulse.



$$\Gamma_{\beta} = n(4)A_{4 \rightarrow 2}$$

Γ_{β} photons/s/cm³ Population density of state 4 Einstein coefficient for the 4 \rightarrow 2 radiative decay rate



$$\Gamma_{\beta} = n_H n_e C_{H,\beta}(T_e, n_e) + n_{H_2} n_e C_{H_2,\beta}(T_e, n_e)$$

Hβ Emission rate

Constants determined by collisional radiative model

$$\frac{\Gamma_{\gamma}}{\Gamma_{\beta}} = f\left(T_e, n_e, \frac{n_{H_2}}{n_H}\right)$$

- Collisional Radiative Model:

$$\begin{aligned} \frac{dn(i)}{dt} = & \left(\overset{\text{Collisional excitation/deexcitation into state } i}{n_e \sum_{k \neq i} n(k) \langle \sigma v \rangle_{k \rightarrow i}} + \overset{\text{Radiative decay into state } i}{\sum_{k > i} A_{k \rightarrow i} n(k)} \right) - \left(\overset{\text{Collisional excitation/deexcitation out of state } i}{n_e \sum_{i \neq k} n(i) \langle \sigma v \rangle_{i \rightarrow k}} + \overset{\text{Radiative decay out of state } i}{\sum_{k < i} A_{i \rightarrow k} n(i)} \right) \\ & - \overset{\text{Ionization of state } i}{n_e n(i) \langle \sigma v \rangle_{i \rightarrow \text{ion}}} + \overset{\text{Ground state H excitation to state } i}{n_H n_e \langle \sigma v \rangle_{H \rightarrow i}} + \overset{\text{Population via H}_2 \text{ dissociation}}{n_{H_2} n_e \langle \sigma v \rangle_{H_2 \rightarrow i}} \end{aligned}$$

---- Input ----


---- Output ----

- Ratio of H & H₂
- Range of Electron Densities
- Range of Electron Temperatures

- Emission Constants for Balmer Series
- Ratio between Emission Rates




H₂ – All Molecular Hydrogen
H – All Atomic Hydrogen




CRM Generated Emission Rates of H gamma and H beta

H₂ – All Molecular Hydrogen
H – All Atomic Hydrogen

Expecting a H_{gamma}-H_{beta} ratio between 0.29 – 0.13




- The graph shows the plasma density over the duration of a PFRC-II pulse.
 - Densification caused by acceleration and heating due to the Rotating Magnetic Field.




- The graphs show electron temperature and the three parameters of the CRM over time around the densification and gas puff stages of a PFRC-II pulse.

Densification




Gas Puff




- The graphs show electron temperature and the three parameters of the CRM over time around the densification and gas puff stages of a PFRC-II pulse.
- Significant issue with the Ratio between Hg/Hb during densification.



Densification



Gas Puff




- Data generated by the monochromator during the PFRC-II pulse.
- H Gamma should happen after H beta since H Gamma since H Gamma is




- The Czerny-Turner Monochromator is in the process of being replaced by a Spectrometer.

- Visible light of various wavelengths will be collected at the same time.
- Eliminates alignment issues with emission measurements.


Monochromator



Spectrometer





- Ocean Optics Spectrometer is now installed and will be collecting data.
- Scripts have been updated to take in data in its new format.
- Data needs to be taken to ensure the spectrometer fixes the alignment issue with emission data.
- Comparing Electron Temperature values from H beta and H alpha ratio to H gamma and H beta ratio.





This work was made possible by funding from the Department of Energy for the Summer Undergraduate Laboratory Internship (SULI) program. This work is supported by the US DOE Contract No. DE-ACo2-09CH11466.





For future testing, the Monochromator will be replaced by an Ocean Optics Spectrometer

